



Audible Blood Glucose Management System



| | |
|--------------------------------------|--|
| NDC # | 8611-501001 |
| Meter Technology | Amperometric- Glucose Oxidase BioSensor |
| Set Up | Fully Audible |
| Test Results | Fully Audible |
| Memory Recall | Fully Audible |
| Blood Glucose Averages | Fully Audible |
| Testing Reminder | Fully Audible |
| Insufficient Blood Warning | Fully Audible |
| Control Test Result | Fully Audible |
| Memory Delete Confirmation | Fully Audible |
| Low Battery Warning | Fully Audible |
| Control Test Indicator | Fully Audible |
| Error Messages | Fully Audible |
| Blood Glucose Measuring Range | 20~600 mg/dL (1.1-33.3 mmol/L) |
| Test Time | 6 seconds |
| Sample Size | .7 micro liters |
| Languages | English & Spanish |
| English/Spanish Selection | External Button |
| Volume Adjustment | Tactile External Buttons with 8 Settings |
| Code Free | No Coding Required |
| Average Results Calculator | 7,14,28,60,90 |
| Meal Indicators | Before & After Meals |
| Time Format | A.M. & P.M. |
| Memory Capacity | 500 Test Results with Date and Time |
| Battery Type | 2 x AAA Batteries |
| Data Transfer | Standard Mini-USB port |
| Data Management Software | Solo V2 DMS |
| Meter Dimensions | 55(W) x 100(L) x 18(H) cm |
| Display Dimensions | 5(L) x 4(W) cm |
| Weight | 75 grams |
| Operating Range | 50-104 F° (10-40° C) |
| Auto Power Off with Strip | 5 minutes |
| Auto Power Off without Strip | 20 seconds |
| Alternate Site Testing (AST) | Yes |
| Strip Vial | Drying Decantent |
| Strip Vial Size | Tube = 30mm; Cap 39mm x 55 mm |
| Control Solution | Hi/Low |
| Test Strip Technology | 4 Leads |
| Warranty | 5 Years |
| Date & Time | Continuous Display |
| Blood Sample Technology | Capillary Action |
| Battery Life | 3000 Tests |
| Storage Case | Compact, Durable, Protective, Case |
| Exterior Finish | Anti-Slip with Integrated Rubber Grommets |